dIFU

Digital Instruction for Use

Nasopharyngeal Airway with Fixed Flange Instructions for use	en
Nasopharyngealtubus mit festem Haltering Gebrauchsanweisung	de
Sonde naso-pharyngée avec collerette fixe Consignes d'utilisation	fr
Via aérea nasofaringea con pestaña fija Instrucciones de uso	es
Tubo nasofaringeo con flangia fissa Istruzioni per l'uso	it
Tubo nasofaringeo com flange fixa Instruções de utilização	pt
Nasofaryngeal luftvägstub med fixerad fläns Bruksanvisning	sv
Nasofaryngeale tube met vaste flens Gebruiksaanwijzing	nl
Nasofaryngeal tube med fast flange Brugervejledning	da
Nenänieluputki, jossa kiinteä pää Käyttöohjeet	fi
Nasofaryngal luftvei med fastsatt flense Pakningsvedlegg	no
Sabit flanlı Nazofaringeal Havayolu Kullanım talimatı	tr
Rurka nosowo-gardłowa z nieruchomym kolnierzem Instrukcja obsługi	pl
Ρινοφαρυγγική αεραγωγία με σταθερή στεφάνη Οδηγίες χρήσης	el
Назофарингеальные дыхательные трубки с неподвижным фланцем Инструкция по применению	ru
Nazofaryngeální dýchací kanyla s pevnou manžetou Návod k použití	cs
固定フランジ付き鼻咽喉頭エアウェイ 取扱説明書	ja
帶固定法蘭的鼻咽喉通氣導管 使用說明	zh
고정 플랜지형 비인두용 에어웨이 사용 설명서	ko

Teleflex®

3 Instructions for use	en
5 Gebrauchsanweisung	de
7 Consignes d'utilisation	fr
9 Instrucciones de uso	es
11 Istruzioni per l'uso	it
13 Instruções de utilização	pt
15 Bruksanvisning	sv
17 Gebruiksaanwijzing	nl
19 Brugervejledning	da
21 Käyttöohjeet	fi
23 Pakningsvedlegg	no
25 Kullanım talimatı	tr
27 Instrukcja obsługi	pl
29 Οδηγίες χρήσης	el
31 Инструкция по применению	ru
33 Návod k použití	cs
35 取扱説明書	ja
37 使用說明	zh
39 사용 설명서	ko

NASOPHARYNGEAL AIRWAY WITH FIXED FLANGE en

PLEASE READ ALL INSTRUCTIONS CAREFULLY BEFORE USE!

Product Description

Fixed Flange Nasopharyngeal Airways consist of a fixed flange and a tube.

The following Fixed Flange Nasopharyngeal Airways are not made with natural rubber latex:

- Robertazzi Nasopharyngeal Airways (green, Mediprene material) – Product numbers beginning with 1231 and 185410
- Robertazzi Nasopharyngeal Airway Kits (green, Mediprene material) – Product numbers 1231 and 1232
- PVC Nasopharyngeal Airways (translucent, clear material) – Product number beginning with 1233/125410
- PVC Nasopharyngeal Airway Kits (translucent, clear material) – Product numbers 1233 and 1234

The following Fixed Flange Nasopharyngeal Airways contain natural rubber latex and should not be used in patient with natural rubber latex allergies!

- Robertazzi Nasopharyngeal Airways (yellow, latex material) – Product number beginning with 1230/185300

Intended Use

To provide a patent airway for a patient with a blockage or potential blockage of the nasal cavity or nasopharynx.

Contraindications

Latex allergy (for product contains natural rubber latex only), traumatic brain injury or central facial fractures, basilar skull or cribriform fractures, anticoagulant therapy, coagulation disorders, transphenoidal surgery, and new postoperative rhinoplasty or septoplasty because it can cause tissue trauma or damage the newly altered structural integrity of the surgical site.

Direction for Use

1) Lubricate the airway with a small amount of water soluble lubricant. Gently insert the airway into nostril until the fixed flange rests against the nostril. Do not force. Gently rotate the airway to ease entry. If insertion is not possible, apply the tube via the contralateral nostril.

2) If a suction catheter is used through the central lumen of the nasopharyngeal airway, ensure correct size and lubricate the outer surface of the suction catheter prior to inserting it into the airway for mucus clearance/removing body fluids purpose. Grasp the flange with one hand while inserting the suction tube with the other hand. **Caution:** Avoid using strong force when introducing the suction catheter into the airway.

3) If face mask is used in conjunction with the nasopharyngeal airway, ensure fitting of face mask does not press against or push the airway into the nostril.

Warning

Some products contain natural rubber latex which may cause allergic reactions. Please refer product label.

CAUTION:

- Do not force the system through the nose, should significant resistance exist.
- Do not attempt insertion where lesions or infection in the nose is observed.

Re-processing of medical devices intended for single use only may result in degraded performance or a loss of functionality. Re-use of single use only medical devices may result in exposure to viral, bacterial, fungal, or prionic pathogens. Validated cleaning and sterilization methods and instructions for reprocessing to original specifications are not available for these medical devices. This product is not designed to be cleaned, disinfected, or sterilized.

Storage Instruction

Keep away from sunlight and keep dry. Do not use if the product sterilisation barrier or its packaging is compromised.

Caution: Federal Law restricts this device to sale by or on the order of a physician. (For USA and Canada only)

Explanation of important symbols and markings on the product label.

SIZE	Product size
Ch.	1 Ch. = 1/3mm
Fr.	1 Fr. = 1/3mm
I.D.	Inner diameter
O.D.	Outer diameter
Qty	Quantity, Number of items
	Not made with natural rubber latex
	This product contains natural rubber latex
	Do not use if package is damaged
	*"Green Dot" Package recycling

Sterilization method: see product label

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AP Services A/S, Aulumvej 7