

# NORSE RESCUE® Cohesive Bandage

Latex Free

---



---

Item No:	CL-001119
Manufacturer:	EMSRUN
Certificates:	CE
Size:	7,5 x 450 cm
Color:	White
Country of origin:	Denmark/China
NSN no:	N/A
Tariff:	30059010
Weight:	20 g

NORSE RESCUE® Cohesive Bandage is an airy bandage with a double adhesive effect, so the bandage sticks to itself and not to the skin. This ensures a fast and easy applying on the patient with no need for closure hooks.

The NORSE RESCUE® Cohesive Bandage is for fixation bandaging, retention of wound dressings and light compression support.