



Universal Transfer Mattress

UNI Transfer Mattress is designed in co-operation with professional paramedic users focusing on ergonomics, patient handling, patient safety and comfort for patient.



UNI Transfer Mattress has been tested by trauma teams to secure efficient handling plus quick, safe and easy transfer of patients including heavy patients.

UNI Transfer Mattress approved according to EN 1865 and EN 1789
Art no: 271 030 001

Safety

Easy to adjust safety belts according to standard EN 1865 and EN 1789

Comfort

Extra comfort by using mix of Sensus memory foam 15 mm and comfort 40 mm V-foam which support patients body weight and deliver a safe and comfort feeling.

Ergonomics

Strong handles to support lift, carry or slide of patient on the mattress using 10 strong handles. It is possible to lower the head for upper airway management. Loading and carry capacity 250 kg.

Suitable to be used in combination with Step by Step system.



Designed to fit on top of most common stretcher designs.

- **Hygienic** full welded easy to clean with most common disinfectants
- **Sliding material / reducing friction** material under the mattress, easy to slide patient over to e.g. hospital bed.
- *Fits together with Germa Step by Step sliding kit.



Accessories: Step by Step
Art no: 15 000 000 117



Specification— Universal Transfer Comfort Mattress

PU- Polyester fibre flame retardant
PP-plate and grass fibre material to reduce friction
Foam Sensus — V-shaped (material that deliver extra comfort and side support to the

Dimensions and weight

Length: 185 x 50 x 6 cm
Weight: 4,7 kg
Handles: 10

Ordering Specification

UNI Transfer Mattress

Art. No.

271030001

